

ENG

Instructions for use:

1,25(OH)₂ VITAMIN D TOTAL ELISA

Catalogue number:

RIS021R

For research use only!

 **BioVendor**
R&D®

BioVendor – Laboratorní medicína a.s.

Karásek 1767/1, 621 00 Brno, Czech Republic

+420 549 124 185

info@biovendor.com

sales@biovendor.com

www.biovendor.com

| | | |
|----|------------------------------------|----|
| 1 | INTENDED USE | 3 |
| 2 | STORAGE, EXPIRATION | 3 |
| 3 | INTRODUCTION | 4 |
| 4 | TEST PRINCIPLE | 5 |
| 5 | PRECAUTIONS | 5 |
| 6 | REAGENT SUPPLIED | 6 |
| 7 | MATERIAL REQUIRED BUT NOT SUPPLIED | 7 |
| 8 | PREPARATION OF REAGENTS | 8 |
| 9 | PREPARATION OF SAMPLES | 10 |
| 10 | ASSAY PROCEDURE | 10 |
| 11 | CALCULATIONS | 12 |
| 12 | PERFORMANCE CHARACTERISTICS | 13 |
| 13 | QUALITY CONTROL | 16 |
| 14 | REFERENCE INTERVALS | 16 |
| 15 | ASSAY PROCEDURE - SUMMARY | 17 |
| 16 | EXPLANATION OF SYMBOLS | 18 |

HISTORY OF CHANGES

| Previous version | Current Version |
|-----------------------------|-----------------|
| ENG.007.A | ENG.008.A |
| "History of changes" added. | |
| Chapter 8.5: New table | |

1 INTENDED USE

Immunoenzymetric assay for the quantitative measurement of 1,25(OH)₂ Vitamin D in serum.
For Research Use Only.

2 STORAGE, EXPIRATION

Before opening or reconstitution, all kits components are stable until the expiry date, indicated on the label, if kept at 2 to 8°C; **except the cartridges which must be stored at room temperature (18°C to 25°C).**

After their reconstitution, **calibrators** are stable for 4 weeks at 2-8°C. For longer storage periods, aliquots should be made and kept at -20°C for 4 months maximum. Avoid subsequent freeze-thaw cycles.

After their reconstitution, **controls** are stable for 3 days at 2-8°C. For longer storage periods, aliquots should be made and kept at -20°C for 1 month maximum. Avoid subsequent freeze-thaw cycles.

Freshly prepared Working Wash solution should be used on the same day.

After its use, discard working HRP conjugate.

Use freshly prepared extraction solvent and washing solvent, do not store them.

Alterations in physical appearance of kit reagents may indicate instability or deterioration.

3 INTRODUCTION

3.1 Biological activities

Vitamin D is mainly synthesized in the skin from 7-dehydrocholesterol and is partially from dietary and supplementation origin. In the liver, Vitamin D is hydroxylated on carbon 25 to produce the intermediate 25OH Vitamin D. 25OH Vitamin D is further metabolized before it can carry out the functions of Vitamin D on intestine, kidneys, bone and other organs and tissues. This subsequent reaction takes place in the kidneys and in other tissues. Thus 25OH Vitamin D is further hydroxylated in the 1α -position to produce $1,25(\text{OH})_2$ Vitamin D. In addition to the above-mentioned tissues, placenta of pregnant women and macrophage cells in case of sarcoidis can also produce some amount of $1,25(\text{OH})_2$ Vitamin D.

$1,25(\text{OH})_2$ Vitamin D is the active form of Vitamin D with regard to the known functions whereas 25OH Vitamin D and Vitamin D itself can be excluded as being physiologically functional. $1,25(\text{OH})_2$ Vitamin D stimulates the intestinal absorption of both calcium and phosphorus. It also stimulates bone resorption and mineralization thereby preventing the development of rickets and osteomalacia.

$1,25(\text{OH})_2$ Vitamin D is also be active in other tissues responsible for Calcium transport (placenta, kidney, mammary gland,...) and endocrine glands such as parathyroid glands. $1,25(\text{OH})_2$ Vitamin D is rapidly metabolized and its half-life is approximately 12h in plasma. Its main metabolite is calcitric acid, a C-23 carboxylic derivative, essentially without any biological activity. In addition to this pathway, $1,25(\text{OH})_2$ Vitamin D undergoes 24-hydroxylation to produce $1,24,25$ -trihydroxyvitamin D. This compound has less biological activity than its parent and this metabolic route is considered as a minor pathway.

The levels of $1,25(\text{OH})_2$ Vitamin D in plasma or serum is 100 to 1000 less than that of 25OH Vitamin D.

Due to its low concentrations and the presence of many similar metabolites, the measurement of $1,25(\text{OH})_2$ Vitamin D requires extraction and separation by chromatography.

3.2 Clinical application

The measurement of circulating $1,25(\text{OH})_2$ Vitamin D is indicated in several disorders affecting calcium metabolism such as: phosphate diabetes, sarcoidosis, renal failure, hyper and hypoparathyroidism, rickets, tumor-associated hypercalcemia, hypercalciuria, Vitamin-resistant dysfunction and treatment with anticonvulsive medication.

4 TEST PRINCIPLE

Only samples and controls, not the calibrators, are extracted with a mixture of solvents and applied on cartridges to separate 1,25(OH)₂ Vitamin D from the other Vitamin D metabolites. After elution of the 1,25(OH)₂ Vitamin D from the samples and controls cartridges, the calibrators, eluted samples and eluted controls are incubated directly in microtiterplate coated with anti-1,25(OH)₂ Vitamin D antibodies.

After an overnight incubation at 4°C, the microtiter plate is washed and the working conjugate solution is added and incubated for 1 hour at 4°C.

The microtiterplate is then washed to stop the competition reaction. The Chromogenic solution (TMB) is added and incubated for 15 minutes at room temperature (18°C to 25°C). The reaction is stopped with the addition of Stop Solution and the microtiterplate is read at the appropriate wavelength.

The amount of 1,25(OH)₂ Vitamin D is determined colourimetrically by measuring the absorbance, which is inversely proportional to the 1,25(OH)₂ Vitamin D concentration.

A calibration curve is plotted and the 1,25(OH)₂ Vitamin D concentrations of the samples are determined by dose interpolation from the calibration curve.

5 PRECAUTIONS

Safety

For research use only

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HBsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections.

Therefore, handling of reagents, serum or plasma specimens should be in accordance with local safety procedures.

All animal products and derivatives have been collected from healthy animals.

Bovine components originate from countries where BSE has not been reported.

Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with all reagents, Stop Solution contains HCl. In case of contact, wash thoroughly with water.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

6 REAGENT SUPPLIED

| Reagents | 96 tests Kit | Color Code | Reconstitution |
|---|------------------------|------------|--|
| Microtiterplate with 96 breakable wells coated with anti-1,25(OH) ₂ Vitamin D antibodies | 96 wells | blue | Ready for use |
| Concentrated HRP | 1 vial 0.2 ml | yellow | Dilute 200 x with conjugate buffer |
| Conjugate Buffer with casein and proclin | 1 vial 30 ml | red | Ready for use |
| 1,25(OH) ₂ Vitamin D Concentrated Conjugate | 1 vial 1 ml | blue | Dilute 40 x with conjugate buffer |
| Incubation Buffer with proclin | 1 vial 20 ml | green | Ready for use |
| Calibrators - N = 0 to 5 (see exact values on vial labels) in phosphate buffer with bovine casein and gentamycin | 6 vials lyophilized | yellow | Add 1.0 ml distilled water |
| Controls - N = 1 or 2 in human plasma with gentamycin | 2 vials lyophilized | silver | Add 3.0 ml distilled water |
| Wash Solution (Tris-HCl) | 1 vial 10 ml | brown | Dilute 200 x with distilled water (use a magnetic stirrer). |
| Chromogen TMB (Tetramethylbenzidine) | 1 vial 25 ml | brown | Ready for use |
| Stop Solution: HCl 1.5N | 1 vial 12 ml | | Ready for use |
| Elution Solution: contains methanol | 1 vial 30 ml | white | Ready for use |
| Adhesive Strips | 4 | | |
| Extraction cartridges | 42 | | Store at R.T. |

Note:

Use Calibrator 0 for dilution of samples with values above the highest calibrator (dilute before extraction step).

7 MATERIAL REQUIRED BUT NOT SUPPLIED

1. Distilled water
2. Diisopropylether ("for analysis"; GC purity \geq 99%)
3. Cyclohexane ("for analysis"; GC purity \geq 99.5 %)
4. Ethyl acetate ("for analysis"; GC purity \geq 99.5 %)
5. Ethanol absolute ("for analysis"; GC purity \geq 99.9 %)
6. Dichloromethane ("for analysis"; GC purity \geq 99.8 %)
7. NB: A Biovendor extraction kit containing all these solvents is available under the reference #RIS024R. This kit contains quantities of solvents necessary to extract the controls and samples, for 2 kits of 1,25(OH)₂ Vitamin D ELISA, in duplicate measurements.
8. Pipettes for delivery of: 50 μ l, 100 μ l, 150 μ l, 200 μ l, 1ml and 2 ml (the use of accurate pipettes with disposable plastic tips is recommended)
9. Glass tubes (12 x 75 mm) for extraction and for elution (closed with a cap for the extraction step).
10. Glass tubes (16 x 100 mm) or (12 x 120 mm), or polypropylene tubes (falcon 2097), for the washing of the cartridges.
11. Vortex mixer
12. Magnetic stirrer
13. Centrifuge operating at 800 g
14. Microtiterplate reader capable of reading at 450 nm and 650 (bichromatic reading)

8 PREPARATION OF REAGENTS

Calibrators

Reconstitute the calibrators with 1.0 ml distilled water.

8.1 Controls

Reconstitute the controls with 3.0 ml distilled water, **carefully to avoid overflow**.

8.2 Working HRP conjugate solution

! The working HRP conjugate solution is to be prepared during the incubation and minimum 1h before its use (cf X.III).

Prepare an adequate volume of working HRP conjugate solution by mixing the 3 reagents in the following sequence: (1) Conjugate buffer, (2) Concentrated Conjugate, (3) Vortex, (4) Concentrated HRP, (5) Vortex.

The order of addition of those 3 reagents is critical and should be rigorously respected to get reproducible Optical Densities.

Prepare the solution according to the number of used strips, as indicated in the below table: for example for 6 strips (48 wells): 250 μ l of concentrated conjugate and 50 μ l of concentrated HRP to 10 ml of conjugate buffer.

Use a vortex to homogenize.

Until its use, keep the working HRP conjugate at room temperature and avoid direct sunlight or use a brown glass vial for its preparation.

The preparation of working HRP conjugate is not stable and must be discarded if not used.

| Number of strips | Volume of Concentrated Conjugate (μ l) | Volume of Concentrated HRP (μ l) | Volume of Conjugate Buffer (ml) |
|------------------|---|---------------------------------------|---------------------------------|
| 1 | 75 | 15 | 3 |
| 2 | 125 | 25 | 5 |
| 3 | 150 | 30 | 6 |
| 4 | 200 | 40 | 8 |
| 5 | 225 | 45 | 9 |
| 6 | 250 | 50 | 10 |
| 7 | 300 | 60 | 12 |
| 8 | 350 | 70 | 14 |
| 9 | 400 | 80 | 16 |
| 10 | 450 | 90 | 18 |
| 11 | 500 | 100 | 20 |
| 12 | 550 | 110 | 22 |

8.3 Working Wash solution

Prepare an adequate volume of Working Wash solution by adding 199 volumes of distilled water to 1 volume of Wash Solution (200x). Use a magnetic stirrer to homogenize. Discard unused Working Wash solution at the end of the day.

8.4 Extraction solvent

2 ml for each control or sample to be tested are needed.

Prepare a fresh solution of diisopropylether, cyclohexane and ethyl acetate: 50/40/10 (volume/volume) according to the number of extractions, as indicated in the table below.

Be careful: the exact proportion of each solvents has to be strictly respected.

| Nb of extraction* | Diisopropylether (ml) | Cyclohexane (ml) | Ethyl acetate (ml) |
|-------------------|-----------------------|------------------|--------------------|
| 1 | 1.1 | 0.9 | 0.2 |
| 8 | 9.2 | 7.4 | 1.8 |
| 16 | 18.4 | 14.7 | 3.7 |
| 42 | 48.3 | 38.6 | 9.4 |

*Patient samples and controls

8.5 Washing solvent:

1 ml for each control or sample to be tested is needed.

Prepare a fresh solution of diisopropylether, cyclohexane, ethyl acetate and ethanol absolute (50/40/10/1 volume/volume) according to the number of extractions, as indicated in the table below.

Be careful: the exact proportion of each solvents has to be strictly respected.

| Nb of extraction* | Diisopropylether (ml) | Cyclohexane (ml) | Ethyl acetate (ml) | Ethanol (µl) |
|-------------------|-----------------------|------------------|--------------------|--------------|
| 1 | 0.6 | 0.5 | 0.1 | 11 |
| 8 | 4.6 | 3.7 | 0.9 | 92 |
| 16 | 9.2 | 7.4 | 1.8 | 184 |
| 42 | 24.1 | 19.3 | 4.8 | 483 |

9 PREPARATION OF SAMPLES

The kit is suitable for serum samples.

Serum samples must be kept at 2-8°C.

If the test is not run within 24 hrs, storage in aliquots, at -20°C is recommended.

Avoid subsequent freeze-thaw cycles.

After thawing, the samples should be vortexed and centrifuged.

10 ASSAY PROCEDURE

10.1 Extraction step: ! Only for controls and samples!

1. Label glass tubes (12x75 mm) for extraction: 2 controls and up to 40 samples.
2. Add 0.5 ml control or sample in the respective tubes.
3. Dispense 2ml extraction solvent in each tube.
4. Tubes are closed with a cap and placed on a shaker for 1 hour at 1200 rpm.
5. Centrifuge each tube for 5 minutes at room temperature (18°C - 25°C at 800 g).
6. Supernatants are needed for the next step of separation.

10.2 Separation step: ! Only for controls and samples!

1. Label glass tubes (16 x 100 mm) or (12 x 120 mm), or polypropylene tubes (falcon 2097), for washing cartridges: 2 controls and up to 40 samples.
2. Put one silica cartridge in each tube.
3. Apply 1.6ml of supernatant (2 x 0.8 ml), obtained after extraction step, on cartridge. Let draw by gravity.
4. Wash cartridges with 1ml washing solvent (cf: reagent preparation).
5. ! Be careful: never apply vacuum on cartridges, just let solvent draw by gravity.
6. Add 500µl dichloromethane on each cartridge, let draw by gravity.
7. Add 500µl of distilled water on each cartridge and centrifuge each tube for 5 minutes at room temperature (18°C - 25°C at 800 g).
8. Label glass tubes (12 x 75 mm) for elution of 1,25(OH)₂ Vitamin D. After centrifugation, transfer cartridges in the corresponding glass tubes.
9. Apply 300µl elution solution on each cartridge to elute 1,25(OH)₂ Vitamin D and centrifuge for 5 minutes at room temperature (18°C - 25°C at 800 g).
10. **Vortex** the eluted fraction.

Note : After this step, samples must be incubated in coated microtiterplate immediately to avoid degradation.

10.3 Incubation step:

1. Select the required number of strips for the run. The unused strips should be resealed in the bag with a desiccant and stored at 2-8°C.
2. Secure the strips into the holding frame.
3. Vortex briefly reconstituted calibrators, extracted controls and extracted samples.
4. Pipette 150 µl of Incubation Buffer into all wells.
5. Pipette 50 µl of each Calibrator (not extracted), eluted controls and eluted samples into the appropriate wells.
6. Incubate for 18±2 hours, at 2-8°C. Cover the plate with a lid or a sealing film.

Prepare the Working HRP conjugate solution 60 min +/- 15 min before washing the wells after the overnight incubation.

1. Aspirate the liquid from each well.
2. Wash the plate 3 times by:
 - dispensing 0.35 ml of Wash Solution into each well and
 - aspirating the content of each well
3. Pipette 200 µl of Working HRP conjugate solution into each well.
4. Incubate for 1 hour at 4°C. Cover the plate with a lid or a sealing film.
5. Aspirate the liquid from each well.
 - Wash the plate 3 times by:
 - dispensing 0.35 ml of Wash Solution into each well and
 - aspirating the content of each well
6. Pipette 200 µl of the Chromogenic solution into each well within 15 minutes following the washing step.
7. Incubate the microtiterplate for 15 minutes at room temperature (18°C to 25°C), avoid direct sunlight.
8. Pipette 100 µl Stop Solution into each well.
9. Read absorbances at 450 nm (reference filter 630 nm or 650 nm) within 1 hour and calculate the results as described in section 12 (CALCULATIONS).

11 CALCULATIONS

Read the plate at 450 nm against a reference filter set at 650 nm (or 630 nm).

Calculate the mean of duplicate determinations.

We recommend the use of computer assisted methods to construct the calibration curve.

4-parameter logistic function curve fitting is the preferred method. Reject obvious outliers.

By interpolation of the sample OD values, determine the 1,25(OH)₂ Vitamin D concentrations of the samples from the calibration curve.

Typical data

The following data are for illustration only and should never be used instead of the real time calibration curve.

| 1,25(OH) ₂ Vitamin D ELISA | OD units |
|---------------------------------------|----------|
| Calibrator: | |
| 0 pg/ml | 2.93 |
| 3 pg/ml | 2.52 |
| 12 pg/ml | 1.85 |
| 50 pg/ml | 1.11 |
| 120 pg/ml | 0.57 |
| 180 pg/ml | 0.36 |

Note: 1 pg/ml = 2.4 pmol/l

Example Version

12 PERFORMANCE CHARACTERISTICS

12.1 Detection Limit

Twenty zero calibrators were assayed along with a set of other calibrators. The detection limit, defined as the apparent concentration two standard deviations below the average OD at zero binding, was 0.8 pg/mL.

12.2 Specificity

Cross reactivity of the 1,25(OH)₂ Vitamin D ELISA assay was determined by testing sera with spiked and unspiked cross reactants. The results are summarized in the following table:

| Compound and Concentration | % Cross reaction |
|---|------------------|
| 1,25(OH) ₂ -Vitamin.D3 at 200 pg/ml | 114 |
| 1,25(OH) ₂ -Vitamin.D2 at 200 pg/ml | 108 |
| 25OH-Vitamin-D3 at 1µg /ml | 0.004 |
| 25OH-Vitamin-D2 at 1µg /ml | 0.0003 |
| 24,25(OH) ₂ -Vitamin.D3 at 200 ng/ml | 0.03 |
| 25,26(OH) ₂ -Vitamin.D3 at 400ng/ml | 0.02 |

The effect of potential interfering substances on samples using the BioVendor 1,25(OH)₂ Vitamin D ELISA test was evaluated. Different levels of Hemoglobin, Bilirubin (conjugated and unconjugated), Triglyceride and Vitamin C in serum samples were tested on samples with different 1,25(OH)₂ Vitamin D Concentration. Our acceptance criteria was to have interference of less than 15%. The tested substances did not affect the performance of the BioVendor 1,25(OH)₂ Vitamin D ELISA.

| Substance | 1,25(OH) ₂ Vitamin D (ng/ml) | Concentration of Interferent (mg/dl) | Mean % Variation |
|------------------------|---|--------------------------------------|------------------|
| Hemoglobin | 31.8 | 250 | 5.0% |
| | | 500 | |
| | 186.5 | 250 | |
| | | 500 | |
| Bilirubin Conjugated | 31.8 | 50 | -12.3% |
| | 186.5 | 50 | |
| Bilirubin Unconjugated | 31.8 | 50 | -0.4% |
| | | 100 | |
| | 186.5 | 50 | |
| | | 100 | |
| Triglyceride | 31.8 | 50 | -1.0% |
| | | 100 | |
| | | 250 | |
| | 186.5 | 50 | |
| | | 100 | |
| | | 250 | |
| Vitamin C | 31.8 | 100 | 4.9% |
| | | 1000 | |
| | 186.5 | 100 | |
| | | 1000 | |

12.3 Precision

| INTRA ASSAY | | | | INTER ASSAY | | | |
|-------------|----|------------------|------|-------------|---|------------------|--------|
| Sample | N | <X> ± SD (ng/ml) | CV % | Sample | N | <X> ± SD (ng/ml) | CV (%) |
| A | 13 | 18.3±2.5 | 13.9 | A | 8 | 26.7±3.5 | 13.2 |
| B | 13 | 168.9±8.4 | 5.0 | B | 8 | 83.4±14.6 | 17.5 |

SD: Standard Deviation; CV: Coefficient of variation

12.4 Accuracy

The sample was diluted with Calibrator 0, before extraction step.

| DILUTION TEST | | | | | | |
|-----------------|-----------------------------------|--------------------------------|-------|-------------|------|--------------|
| Sample Dilution | Theoretical concentration (ng/ml) | Measured concentration (ng/ml) | Slope | Y-Intercept | R | Recovery (%) |
| 1/1 | 118.9 | 118.9 | 0.99 | 1.12 | 0.99 | 100 |
| 1/2 | 59.4 | 60.7 | | | | 102 |
| 1/4 | 29.7 | 29.3 | | | | 99 |
| 1/8 | 14.8 | 16.6 | | | | 112 |

Conversion factor:

From pg/ml to pmol/l: x 2.4

From pmol/l to pg/ml: x 0.42

To the best of our knowledge, no international reference material exists for this parameter.

| RECOVERY TEST | | |
|---|--|---------------|
| Added 1,25(OH) ₂ - Vitamin D (pg/ml) | Recovered 1,25(OH) ₂ -Vitamin D (pg/ml) | Recovered (%) |
| 52.4 | 54.1 | 103 |
| 104.7 | 111.1 | 106 |
| 157.1 | 155.8 | 99 |

13 QUALITY CONTROL

- If the results obtained for Control 1 and/or Control 2 are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots. Controls which contain azide will interfere with the enzymatic reaction and cannot be used.
- Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practises.
- It is recommended that Controls be routinely assayed as unknown samples to measure assay variability. The performance of the assay should be monitored with quality control charts of the controls.
- It is good practise to check visually the curve fit selected by the computer.

14 REFERENCE INTERVALS










These values are given only for guidance; each laboratory should establish its own normal range of values.

Normal samples tested with BioVendor Elisa assay were measured between 19.3 and 53.8 pg/ml. Patients with renal failure (n = 20) were measured < 6,9 pg/ml.

15 ASSAY PROCEDURE - SUMMARY

| | CALIBRATORS μl | SAMPLE(S) CONTROLS μl |
|----------------------------------|--|--------------------------|
| EXTRACTION | | |
| Calibrators | - | - |
| Samples / Controls | - | 500 |
| Extraction solvent | - | 2000 |
| Shaking | 1 hour at 1200 rpm | |
| Centrifugation | 5 minutes at 800 g | |
| SEPARATION | | |
| Supernatant from extraction step | - | 1600 |
| CARTRIDGE | | |
| Supernatant | 1600μl | |
| Washing Solvent | 1000μl | |
| Dichloromethane | 500μl | |
| Distilled water | 500μl | |
| Centrifugation | 5 minutes at 800 g | |
| Elution solution | 300 μl | |
| Centrifugation | 5 minutes at 800 g | |
| | Vortex | |
| INCUBATION STEP | | |
| In microtiterplate | | |
| Incubation Buffer | 150μl | 150μl |
| Calibrators | 50μl | - |
| Extracted samples | - | 50μl |
| | Cover the plate with a lid or sealing film Incubate 18 ± 2 h (overnight) at 4°C (2-8°C) Prepare working HRP solution 1 hour before next step Aspirate the contents of each well Wash 3 times with 350μl of Wash Solution and aspirate | |
| Working HRP Conjugate | 200 μl | 200 μl |
| | Cover the plate with a lid or sealing film and Incubate for 1 hour at 4°C (2-8°C). Aspirate the contents of each well. Wash 3 times with 350 μl of Wash Solution and aspirate. | |
| TMB | 200 μl | 200 μl |
| | Incubate for 15 min at room temperature (18°C to 25°C). | |
| Stop Solution | 100 μl | 100 μl |
| | Read on a microtiterplate reader. Record the absorbance of each well at 450 nm (versus 630 or 650 nm). | |

16 EXPLANATION OF SYMBOLS

| | |
|--|---|
|  | Catalogue number |
|  | Batch code |
|  | Caution |
|  | Use by date |
|  | Temperature limit |
|  | Manufacturer |
|  <p data-bbox="256 1184 464 1216">www.biovendor.com</p> | Read electronic instructions for use - eIFU |
|  | The content is sufficient for 96 tests |
|  | Biological risks |



BioVendor – Laboratorní medicína a.s.
Karásek 1767/1, 621 00 Brno, Czech Republic
+420 549 124 185
info@biovendor.com
sales@biovendor.com
www.biovendor.com

Example Version

