

ENG

Instructions for Use:
CANINE CYSTATIN C ELISA

Catalogue number:
RBL021R

For research use only!

 **BioVendor**
R&D[®]



BioVendor – Laboratorní medicína a.s.

Karásek 1767/1, 621 00 Brno, Czech Republic

+420 549 124 185

info@biovendor.com

sales@biovendor.com

www.biovendor.com

1. INTENDED USE	3
2. STORAGE, EXPIRATION	3
3. INTRODUCTION	3
4. TEST PRINCIPLE	4
5. PRECAUTIONS	4
6. REAGENT SUPPLIED	4
7. MATERIAL REQUIRED BUT NOT SUPPLIED	4
8. PREPARATION OF REAGENTS	5
9. PREPARATION OF SAMPLES	6
10. ASSAY PROCEDURE	6
11. PERFORMANCE CHARACTERISTICS	7
12. DEFINITION OF THE STANDARD	8
13. REFERENCE RANGE	8
14. CALCULATION	9
15. REFERENCES	10
16. EXPLANATION OF SYMBOLS	11
17. ASSAY PROCEDURE - SUMMARY	12

HISTORY OF CHANGES

Previous version	Current version
	ENG.001.A
New edition	

1. INTENDED USE

Enzyme Immunoassay for the quantitative determination Cystatin C in canine serum and canine urine.

2. STORAGE, EXPIRATION

The kit must be stored at 2 – 8°C.

The opened components can be stored for one week at 2 – 8°C.

3. INTRODUCTION

Cystatin C is a non-glycosylated basic protein belonging to the super-family of cysteine proteinase inhibitors. It consists of a single polypeptide chain having 120 amino acids.

It is produced by all nucleated cells within the body and is released during phagocytosis and inflammation. In the kidney, cystatin C is freely filtrated through the glomerulus and reabsorbed and catabolized in the proximal renal tubules. The rate of cystatin C synthesis is constant, independent of age, gender and muscle mass. High concentrations can be found in serum, seminal fluid, cerebrospinal fluid (CSF), and synovial fluid, and lower concentrations can be found in urine.

In human medicine, cystatin C is the most important endogenous serum marker of renal function assessment. Cystatin C evaluation is able to detect an earlier stage of decreased glomerular filtration rate (GFR) than other parameters (serum creatinine, creatinine clearance etc.) and it is considered particularly useful in patients with a high risk of developing nephropathies. Imbalance between cystatin C and cysteine proteinases is associated with inflammation, cancer, Alzheimer's disease, multiple sclerosis and hereditary cystatin C amyloid angiopathy. An increased level has been found in patients with autoimmune diseases. On the other hand, low concentration of cystatin C presents a risk factor for secondary cardiovascular events.

In veterinary medicine, there are multiple reports of the use of cystatin C in the evaluation of renal function indicating that cystatin C is also the most important serum (urine) marker of renal function assessment in dogs.

4. TEST PRINCIPLE

The microtiter plate is coated with the antibody specifically binding the canine Cystatin C. The canine serum or urine is incubated in the plate with the capture antibody.

The specimen is washed out and the specifically bound protein is incubated with HRP-labelled detection antibody. Unbound reagent is then washed out. Horseradish peroxidase (HRP) bound in the complex reacts with the chromogenic substrate (TMB) creating the blue colour. The reaction is stopped by addition of STOP solution (H₂SO₄).

The absorbance values are measured at 450 nm (optionally 450/630 nm) and are proportional to the concentration of Canine Cystatin C in the specimen. The concentration of Cystatin C in unknown samples is determined from the calibration curve which is created by plotting the absorbance values against the standard concentration values.

5. PRECAUTIONS

1. For research use only
2. For professional laboratory use
3. The reagents with different lot numbers should not be mixed
4. To prevent cross sample contamination, use disposable labware and pipette tips
5. To protect laboratory stuff, wear protective gloves and protective clothing
6. The substrate solution should remain colourless, keep it protected from light
7. The test should be performed at standard laboratory conditions (temperature 25°C ±2°C).

6. REAGENT SUPPLIED

Item	Qty.
Antibody Coated Microtiter Plate	96 wells
Conjugate Diluent	13 ml
Ab-HRP Conjugate 50x conc.	0.26 ml
Master Standard (lyophilized)	2 vials
Quality Control A (lyophilized)	2 vials
Quality Control B (lyophilized)	2 vials
Dilution Buffer 5x conc.	13 ml
Wash Buffer 15x conc.	50 ml
Substrate Solution	13 ml
STOP Solution	13 ml

7. MATERIAL REQUIRED BUT NOT SUPPLIED

- Glassware and test tubes
- Microtiter plate washer
- Precision pipettes (various volumes) with tips
- Orbital shaker
- Microtiter plate reader capable of measuring absorbance at 450 nm or 450/630 nm with software for data generation

8. PREPARATION OF REAGENTS

Use new pipette tip for pipetting different reagents and samples to prevent cross-contamination. All reagents and samples should be allowed to reach the temperature 25°C ±2°C.

8.1 Preparation of Standards

Canine Cystatin C Master Standard

Refer to the Certificate of Analysis for current volume of Dilution Buffer needed for reconstitution of standard!!! Reconstitute the lyophilized Master Standard with Dilution Buffer just prior to the assay. Let it dissolve at least 15 minutes with occasional gentle shaking (not to foam). The resulting concentration of the Canine Cystatin C in the stock solution is 10 ng/ml.

Prepare set of Standard solution as follows:

Use the Master Standard to produce a dilution series (as below). Mix each tube thoroughly before the next transfer. The Dilution Buffer serves as Blank.

	Volume of Standard	Dilution Buffer	Concentration
Std1	Standard 25 ng/ml (lyophilised)	See CoA	10 ng/ml
Std2	300 µl of Std1	300 µl	5 ng/ml
Std3	300 µl of Std2	300 µl	2.5 ng/ml
Std4	300 µl of Std3	300 µl	1.25 ng/ml
Std5	300 µl of Std4	300 µl	0.625 ng/ml
Std6	300 µl of Std5	300 µl	0.31 ng/ml
Blank	-	200 µl	0 ng/ml

8.2 Dilution Buffer 1x

Dilute Dilution Buffer Concentrate (5x) five-fold in distilled water to prepare a 1x working solution. Example: 13 ml of Dilution Buffer Concentrate (5x) + 52 ml of distilled water for use of all 96 wells.

Stability and storage:

The diluted Dilution Buffer is stable 1 week when stored at 2-8°C. Opened Dilution Buffer Concentrate (5x) is stable 3 months when stored at 2-8°C.

8.3 Quality Control A and B

Reconstitute the lyophilized Quality Controls in dilution buffer, for the volume information see the Certificate of Analysis. Let the QCs rehydrate for 15 min. The reconstituted Quality controls are ready to use, do not dilute them.

8.4 Wash Buffer 1x

Prepare a working solution of Wash Buffer by adding 50 ml of Wash Buffer 15x conc. to 700 ml of deionized/ distilled water (dH₂O). Mix well. Store at 4°C for two weeks or at -20°C for long term storage.

8.5 Ab-HRP Conjugate 1x

Prepare a working solution of Ab-HRP by adding 265 µl of Ab-HRP Conjugate 50x conc. to 13 ml of Conjugate Diluent. Mix well. The diluted HRP-labelled antibody is stable for two weeks when stored at 2-8°C.

9. PREPARATION OF SAMPLES

Canine serum or urine may be used with this assay. It is recommended to assay not-frozen samples. For long-term storage the serum and plasma samples should be frozen at minimum -70°C.

Recommended dilution of **serum** samples is 1:1000. It is recommended to use the two-step dilution.

Dilution A (20x) for both singlets and duplicates: 5 µl of samples + 95 µl of Dilution Buffer.

Dilution B (50x): 10 µl of Dilution A + 490 µl of Dilution Buffer, for both singlets and duplicates; Mix well, vortex is recommended.

Recommended dilution of **urine** samples is 1:100. Add 5 µl of samples to 495 µl of Dilution Buffer. Dilute urine samples from dogs with serum creatinine > 1.4 mg/dl at least 1:1000. Two step dilution is recommended.

Dilution A (20x) for both singlets and duplicates: 5 µl of samples + 95 µl of Dilution Buffer.

Dilution B (50x): 10 µl of Dilution A + 490 µl of Dilution Buffer, for both singlets and duplicates; Mix well, vortex is recommended.

Do not store the diluted samples. Do not store the diluted samples.

10. ASSAY PROCEDURE

1. Prepare the reagents as described in the previous chapter.
2. Prepare the reagents as described in the previous chapter.
3. Pipette 100 µl of set of Standards, Quality Controls, diluted Samples and Dilution Buffer = Blank into each well. Incubate for **1 hour** at 25°C ±2°C, shaking at 300 rpm.
4. Wash the wells 3-times with 1x Wash Buffer (350 µl/well). When finished, tap the plate against the paper towel to remove the liquid completely.
5. Pipette 100 µl of HRP-labelled Antibody Conjugate into each well. Incubate for **1 hour** at 25°C ±2°C, shaking at 300 rpm.
6. Wash the wells as described in point 3.
7. Pipette 100 µl Substrate solution, incubate for **15 min**, at 25°C ±2°C. Avoid exposure to the light during this step.
8. Pipette 100 µl of STOP solution.
9. Read the signal at 450 or 450/630 nm within 15 min.

11. PERFORMANCE CHARACTERISTICS

Samples used in the tests were diluted 1:1000 as recommended and assayed. The results are multiplied by the dilution factor.

11.1 Sensitivity

The limit of detection, defined as a concentration of canine Cystatin C giving absorbance higher than absorbance of blank + 3 standard deviations, is better than 0,005 ng/ml of sample.

11.2 Precision

Intra-assay (Within-Run) (n=8)

Sample	Mean (µg/ml)	SD (µg/ml)	CV (%)
1	0.68	0.06	8.4
2	4.60	0.19	4.2

Inter-assay (Run-to-run) (n=8)

Sample	Mean (µg/ml)	SD (µg/ml)	CV (%)
1	0.69	0.03	4.1
2	4.42	0.32	7.3

11.3 Accuracy

11.3.1 Dilution linearity

Sample	Dilution	Measured concentration (µg/ml)	Expected concentration (µg/ml)	Yield (%)
Serum 1		2.88	-	-
	2x	1.36	1.44	95
	4x	0.64	0.72	88
	8x	0.31	0.36	86
Serum 2		7.77	-	-
	2x	3.60	0.39	93
	4x	1.71	1.94	88
	8x	0.86	0.97	88
Urine 1		0.74	-	-
	2x	0.36	0.37	96
	4x	0.17	0.18	92
	8x	0.09	0.09	93
Urine 2		12.55	-	-
	2x	5.97	6.28	95
	4x	2.88	3.14	92
	8x	1.42	1.57	91

11.3.2 Spiking Recovery

Sample	Measured concentration (µg/ml)	Expected concentration (µg/ml)	Yield (%)
Serum 1	1.52	-	-
	2.10	2.02	104
	2.62	2.52	104
	3.63	3.52	103
Serum 2	1.05	-	-
	1.66	1.55	108
	2.11	2.05	103
	3.16	3.05	104
Urine 1	0.51	-	-
	0.60	0.57	106
	0.63	0.63	100
	0.79	0.76	104
Urine 2	0.16	-	-
	0.22	0.22	100
	0.30	0.28	107
	0.42	0.41	102

12. DEFINITION OF THE STANDARD


The recombinant canine Cystatin C is used as the Standard. The recombinant canine Cystatin C, produced in E.coli, is 14.85 kDa protein.

13. REFERENCE RANGE

It is recommended that each laboratory include its own panel of control samples in the assay. Each laboratory should establish its own normal and pathological reference ranges for canine Cystatin C levels with the assay.

14. CALCULATION










The standard curve needs to be measured in every test. Most of the microplate reader can automatically calculate the analyte concentration using 4-parameter algorithm or alternative functions to fit the standard points properly. The concentrations need to be multiplied by the dilution factor, either automatically by reader or manually.



15. REFERENCES

1. Monti P, Benchekroun G, Berlato D, Archer J: Initial evaluation of canine urinary cystatin C as a marker of renal tubular function. *J Small Anim Pract*; 53(5):254-259 (2012)
2. Kavitha K, Yathiraj S, Ramachandra SG: Serum cystatin-C as a marker for renal dysfunction and its correlation with creatinine and blood urea nitrogen (BUN). *JCVA*; 27(1):15-17 (2011)
3. Jonkisz P, Kungl K, Sikorska A, Kurosad A, Nicpoń J: Cystatin C analysis in the dog: a comparison of turbidimetric and nephelometric assay results. *Acta Vet Hung*; 58(1):59-67 (2010)
4. Miyagawa Y, Takemura N, Hirose H: Evaluation of the Measurement of Serum Cystatin C by an Enzyme-Linked Immunosorbent Assay for Humans as a Marker of the Glomerular Filtration Rate in Dogs. *J Vet Med Sci*; 71(9):1169-1176 (2009)
5. Pasa S, Kilic N, Atasoy A, Derincegoz OO, Karul A: Serum Cystatin C Concentration as a Marker Acute Renal Dysfunction in Critically Ill Dogs. *J Anim Vet Adv*; 7(11):1410-1412 (2008)
6. Wehner A, Hartmann K, Hirschberger J: Utility of serum cystatin C as a clinical measure of renal function in dogs. *J Am Anim Hosp Assoc*; 44(3):131-138 (2008)
7. Antognoni MT, Siepi D, Porciello F, Rueca F, Fruganti G: Serum Cystatin-C Evaluation in Dogs Affected by Different Diseases Associated or Not with Renal Insufficiency. *Vet Res Commun*; 31(Suppl. 1): 269-271 (2007)
8. Antognoni MT, Siepi D, Porciello F, Fruganti G: Use of Serum Cistatin C Determination as a Marker of Renal Function in the Dog. *Vet Res Commun*; 29(Suppl. 2):265-267 (2005)
9. Filler G, Bokenkamp A, Hofmann W, Le Bricond T, Martinez-Bru C, Grubb A: Cystatin C as a marker of GFR-history, indications, and future research. *Clin Biochem*; 38:1-8 (2005)
10. Mares J, Stejskal D, Vavrouskova J, Urbanek K, Herzig R, Hlustik P: Use of Cystatin C Determination in Clinical Diagnosis. *Biomed Papers*; 147(2):177-180 (2003)
11. Almy FS, Christopher MM, King DP, Brown SA: Evaluation of cystatin C as an endogenous marker of glomerular filtration rate in dogs. *J Vet Intern Med*; 16(1):45-51 (2002)
12. Coll E, Botey A, Alvarez L, Poch E, Quinto L, Saurina A, Vera M, Piera C, Darnell A: Serum Cystatin C as a New Marker for Noninvasive Estimation of Glomerular Filtration Rate and as a Marker for Early Renal Impairment. *Am J Kidney Dis*; 36(1):29-34 (2000)

16. EXPLANATION OF SYMBOLS

	Catalogue number
	Batch code
	Caution
	Use by date
	Temperature limit
	Manufacturer
 <p data-bbox="256 1182 464 1205">www.biovendor.com</p>	Read electronic instructions for use - eIFU
	The content is sufficient for 96 tests
	Biological risks

17. ASSAY PROCEDURE - SUMMARY

Add 100 μ L of Standards, diluted QCs and Samples to the wells



Incubate 1 hour at 25°C, shaking at 300 rpm

3-times wash the wells (350 μ L/well)



Add 100 μ L of HRP-conjugated Antibody to the wells



Incubate for 1 hour at 25°C, shaking at 300 rpm

3-times wash the wells (350 μ L/well)



Add 100 μ L of Substrate Solution to the wells



Incubate for 15 min in the dark at 25°C, NO shaking

Add 100 μ L of Stop Solution to the wells



Read the signal at 450 nm (450/630 nm) within 15 min



BioVendor – Laboratorní medicína a.s.

Karásek 1767/1, 621 00 Brno, Czech Republic

+420 549 124 185

info@biovendor.com

sales@biovendor.com

www.biovendor.com

