

PRODUCT DATA SHEET

Decorin Human, Rabbit Polyclonal Antibody

Cat. No.: RD181324100

Source of Antigen: E. coli

Host: Rabbit

Size: 0.1 mg

Type: Polyclonal Antibody

Other Names:

Bone proteoglycan II, PG-S2, PG40, DCN, SLRR1B

Preparation:

The antibody was raised in rabbits by immunization with the recombinant Human Decorin.

Amino Acid Sequence of Immunogen:

The immunization antigen (37.6 kDa) is a protein containing 339 AA of recombinant Human Decorin. UniProtKB acc.no. P07585 (Asp31-Lys359). N-terminal His-tag, 10 extra AA (highlighted).

MKHHHHHHAS DEASGIGPEV PDDRDFEPSL GPVCPFRCQC HLRVVQCSDL
GLDKVPKDLP PDTLLDLQN NKITEIKDGD FKNLKLNHAL ILVNNKISKV SPGAFTPLVK
LERLYLSKNQ LKELPEKMPK TLQELRAHEN EITKVRKVTF NGLNQMIVIE LGTNPLKSSG
IENGAFQGMK KLSYIRIADT NITSIPQGLP PSLTELHLDG NKISRVDAAS LKGLNNLAKL
GLSFNSISAV DNGSLANTPH LRELHLDNNK LTRVPGGLAE HKYIQVVYLH NNNISVVGSS
DFCPPGHNTK KASYSGVSLF SNPVQYWEIQ PSTFRCVYVR SAIQLGNYK

Purification Method:

Immunoaffinity chromatography on a column with immobilized recombinant Human Decorin

Species Reactivity:

Human

Antibody Content:

0.1 mg (determined by BCA method, BSA was used as a standard)

Formulation:

The antibody is lyophilized in 0.05 M phosphate buffer, 0.1 M NaCl, pH 7.2. ****AZIDE FREE****.

Reconstitution:

Add 0.2 ml of deionized water and let the lyophilized pellet dissolve completely. Slight turbidity may occur after reconstitution, which does not affect activity of the antibody. In this case clarify the solution by centrifugation.

Shipping:

At ambient temperature. Upon receipt, store the product at the temperature recommended below.

Storage/Stability:

The lyophilized antibody remains stable and fully active until the expiry date when stored at -20°C. Aliquot the product after reconstitution to avoid repeated freezing/thawing cycles and store frozen at -80°C. Reconstituted antibody can be stored at 4°C for a limited period of time; it does not show decline in activity after one week at 4°C.

Quality Control:

Indirect ELISA – to determine titer of the antibody

SDS PAGE – to determine purity of the antibody

BCA - to determine quantity of the antibody

Applications:

ELISA

Note:

This product is for research use only.