

PRODUCT DATA SHEET

Adiponectin Globular Human E. coli

Cat. No.: RD172112100

Size: 0.1 mg

Source: E. coli

Type: Recombinant protein

Species: Human

Other Names:

Adipocyte C1q and collagen domain-containing protein, Adipocyte complement-related 30 kDa protein, ACRP30, Adipose most abundant gene transcript 1 protein, apM-1, Gelatin-binding protein, ADIPOQ, ACDC, APM1, GBP28

Description:

Total 171AA. N-terminal His-tag and the protease cleavage site, 31 extra AA. The AA sequence is homologous to the amino acid sequence (105–244) of the Human full-length Adiponectin.

Amino Acid Sequence:

MSWWHHHHHH NWNIPPTQDT TQDLWFEGAM GGEGAYVYRS AFSVGLETYV TIPNMPPIRFT
KIFYNQQNHY DGSTGKFHCN IPGLYYFAYH ITVYMKDVKV SLFKKDKAML FTYDQYQENN
VDQASGSVLL HLEVGDQVWL QVYGEGERNG LYADNDNDST FTGFLLYHDT N

Purity:

> 90 % by SDS-PAGE

Endotoxin:

< 1.0 EU/μg

Formulation:

Filtered (0.4 μm) and lyophilized in 0.5 mg/mL in 0.05 M phosphate buffer, 0.075 M NaCl, pH 7.4

Reconstitution:

Add deionized water to prepare a working stock solution of approximately 0.5 mg/mL and let the lyophilized pellet dissolve completely.

Shipping:

At ambient temperature. Upon receipt, store the product at the temperature recommended below.

Storage/Stability:

Store the lyophilized protein at -80 °C. Lyophilized protein remains stable until the expiry date when stored at -80 °C. Aliquot reconstituted protein to avoid repeated freezing/thawing cycles and store at -80 °C for long term storage. Reconstituted protein can be stored at 4 °C for a week.

Quality Control:

BCA to determine quantity of the protein.
SDS PAGE to determine purity of the protein.
Endotoxin level determination.

Applications:

ELISA, Western blotting

Note:

This product is intended for research use only.