

PRODUCT DATA SHEET

Macrophage Colony Stimulating Factor Mouse E. coli

Cat. No.: RP1764390010

Size: 10 µg

Source: E. coli

Type: Recombinant protein

Species: Mouse

Other Names:

Description:

Macrophage Colony Stimulating Factor Mouse Recombinant produced in E.coli is a disulfide linked homodimer, non-glycosylated, polypeptide chain containing 2×156 amino acids and having a total molecular mass of 36.4 KD. MCSF is purified by proprietary chromatographic techniques.

Amino Acid Sequence:

MKEVSEHCSH MIGNGHLKVL QQLIDSQMET SCQIAFEFVD QEQLDDPVCY LKKAFFLVQD
IIDETMRFKD NTPNANATER LQELSNLNS CFTKDYEEQN KACVRTFHET PLQLLEKIKN
FFNETKNLLE KDWNIFTKNC NNSFAKCSSR DVVTKP

Purity:

Greater than 98.0% as determined by: (a) Analysis by RP-HPLC. (b) Analysis by SDS-PAGE.

Formulation:

The protein (1 mg/ml) was lyophilized with 0.5×PBS, pH 8.0.

Shipping:

At ambient temperature. Upon receipt, store the product at the temperature recommended below.

Storage/Stability:

Lyophilized M-CSF although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution MCSF should be stored at 4°C between 2–7 days and for future use below -18°C. Please prevent freeze-thaw cycles.