

PRODUCT DATA SHEET

SARS-CoV-2 Spike Glycoprotein S1, Insect Cells

Cat. No.: RN0012500

Size: 0.5 mg

Source: Insect Cells

Type: Recombinant protein

Species: SARS-CoV-2

Other Names:

Severe acute respiratory syndrome coronavirus 2 spike glycoprotein S1, 2019 novel coronavirus S1 protein, SARS-CoV-2 S1 subunit, COVID-19

Description:

SARS-CoV-2 Spike S1, His-tag (Wuhan-Hu-1 strain). Contains Spike protein 1-674, C-terminally tagged with His-tag. Recombinant SARS-CoV-2 Spike S1, His-tag is manufactured in insect cells and purified from culture supernatant by Protein G chromatography.

Formulation:

Presented in Dulbecco's phosphate buffered saline (DPBS) pH 7.4.

Reconstitution:

Defrost at ambient temperature

Shipping:

On ice. Upon receipt, store the product at the temperature recommended below.

Storage/Stability:

Store protein at -80°C. Protein remains stable until the expiry date when stored at -80°C. Avoid repeated freezing/thawing cycles.

Applications:

COVID-19, ELISA, Western blotting

Note:

This product is intended for research use only.