

## PRODUCT DATA SHEET

### Chorionic Gonadotropin Human

**Cat. No.:** RP17825025**Size:** 2.5 mg**Source:** Urine of pregnant women**Type:** Native protein**Species:** Human**Other Names:**

hCG

**Description:**

Human Chorionic Gonadotropin is produced from a sterile preparation of placental glucoprotein urine of pregnant women having a total molecular mass of 36,700 Dalton. The hCG consists of 237 amino acids, a chain-92 amino acids and b chain-145 amino acids. The hCG is purified by proprietary chromatographic techniques.

**Formulation:**

The hCG was lyophilized with no additives.