

PRODUCT DATA SHEET

Pregnant Mare Serum Gonadotropin (PMSG)

Cat. No.: RP17827210000**Size:** 100,000 IU**Source:** Serum of pregnant mares**Type:** Native protein**Other Names:**

PMSG

Description:

PMSG is a complex glycoprotein obtained from the serum of pregnant mares. This 43–63 kDa protein is capable of supplementing and being substituted for the follicle stimulating gonadotrophin and interstitial cell-stimulating hormone of the anterior pituitary gland in both the male and female. Thus PMSG-Intervet stimulates development of the ovarian follicle in the female.

Formulation:

The PMSG was lyophilized with no additives.

Shipping:

At ambient temperature. Upon receipt, store the product at the temperature recommended below.

Storage/Stability:

Lyophilized PMSG although stable at room temperature for 3 weeks, should be stored between 2–8°C.