

## Mouse Anti Adenovirus Hexon Protein Antibody, Clone: A15

### Product Data Sheet

**Host:** Mouse  
**Clone:** A15  
**Isotype:** Mouse IgG

**Cat. No.:** RM812137100 (0.1 mg)

### Species Reactivity

This antibody is specific for the hexon protein of Adenovirus. The antibody has been shown to recognise adenovirus from all strains that have been tested, including adenovirus 2, 3, 4, 5, 15, 31, 40 and 41.

The antibody does not cross react with Astrovirus, *Campylobacter spp.*, *Clostridium difficile*, *Cryptosporidium*, Enterovirus, *Escherichia coli*, *Giardia*, *Helicobacter pylori*, HAV, *Influenza type A*, *Influenza type B*, *Listeria monocytogenes*, Norovirus, Respiratory syncytial virus, Rotavirus, *Salmonella spp.*, *Shigella spp.*, *Yersinia enterocolitica*.

### Purification Method

The antibody was purified by affinity chromatography on protein A.

### Antibody Content

0.1 mg; concentration: 1.0 mg/ml

### Formulation

Supplied as a liquid in Phosphate Buffered Saline pH 7.4, 0.09% Sodium Azide.

### Shipping

At ambient temperature.

### Storage/Stability

Store at +4°C for up to three months, or at -20°C for longer. Avoid repeated freeze/thaw cycles.

### Expiration

See vial label.

### Lot Number

See vial label.

### Applications

ELISA

### Note

For research use only.

<b>HEADQUARTERS:</b>				
BioVendor Laboratorní medicína, a.s.	Karasek 1767/1	621 00 Brno CZECH REPUBLIC	Phone: +420-549-124-185 Fax: +420-549-211-460	E-mail: info@biovendor.com sales@biovendor.com Web: www.biovendor.com
AUSTRIA: BioVendor GesmbH	Nußdorfer Straße 20/10	1090 Vienna AUSTRIA	Phone: +43-1-89090-25 Fax: +43-1-89090-2515	E-mail: infoAustria@biovendor.com
GERMANY, SWITZERLAND: BioVendor GmbH	Otto-Hahn-Straße 16	34123 Kassel GERMANY	Phone: +49-6221-433-9100 Fax: +49-6221-433-9111	E-mail: infoEU@biovendor.com
USA, CANADA AND MEXICO: BioVendor LLC	128 Bingham Rd. Suite 1300	Asheville, NC 28806 USA	Phone: +1-828-575-9250 +1-800-404-7807 Fax: +1-828-575-9251	E-mail: infoUSA@biovendor.com