

## Mouse Anti Influenza B Antibody, Clone: YB91

### Product Data Sheet

**Host:** Mouse  
**Clone:** YB91  
**Isotype:** Mouse IgG<sub>2</sub>

**Cat. No.:**  
RM812164100 (0.1 mg)

### Species Reactivity

This antibody is specific for *Influenza B* detection. The antibody does not cross react with Adenovirus or Respiratory Syncytial virus.

### Purification Method

The antibody was purified by affinity chromatography on protein A.

### Antibody Content

0.1 mg; concentration: 1.0 mg/ml

### Formulation

Supplied as a liquid in Phosphate Buffered Saline pH 7.4, 0.09% Sodium Azide.

### Shipping

At ambient temperature.

### Storage/Stability

Store at +4°C for up to three months, or at -20°C for longer. Avoid repeated freeze/thaw cycles.

### Applications

ELISA

### Note

For research use only.

<b>HEADQUARTERS:</b> BioVendor Laboratorní medicína, a.s.	Karasek 1767/1	621 00 Brno CZECH REPUBLIC	Phone: +420-549-124-185 Fax: +420-549-211-460	E-mail: info@biovendor.com sales@biovendor.com Web: www.biovendor.com
AUSTRIA: BioVendor GesmbH	Nußdorfer Straße 20/10	1090 Vienna AUSTRIA	Phone: +43-1-89090-25 Fax: +43-1-89090-2515	E-mail: infoAustria@biovendor.com
GERMANY, SWITZERLAND: BioVendor GmbH	Otto-Hahn-Straße 16	34123 Kassel GERMANY	Phone: +49-6221-433-9100 Fax: +49-6221-433-9111	E-mail: infoEU@biovendor.com
USA, CANADA AND MEXICO: BioVendor LLC	128 Bingham Rd. Suite 1300	Asheville, NC 28806 USA	Phone: +1-828-575-9250 +1-800-404-7807 Fax: +1-828-575-9251	E-mail: infoUSA@biovendor.com