

Mouse Anti Respiratory Syncytial Virus Antibody, Clone: RV11

Product Data Sheet

Host: Mouse
Clone: RV11
Isotype: Mouse IgG₁

Cat. No.: RM812398100 (0.1 mg)

Species Reactivity

This antibody is specific for respiratory syncytial virus. The antibody does not cross react with Influenza Type A, Influenza Type B or Adenovirus.

Purification Method

Purified in a multi-step procedure to ≥95% purity by SDS-PAGE.

Antibody Content

0.1 mg; concentration: 1.0 mg/ml

Formulation

Supplied as a liquid in Phosphate Buffered Saline pH 7.4, 0.09% Sodium Azide.

Shipping

At ambient temperature.

Storage/Stability

Short term storage at 2-8°C.
Longer term storage at -20°C.
Avoid repeated freeze/thaw cycles.

Applications

ELISA, LFA

Note

For research use only.

HEADQUARTERS:				
BioVendor Laboratorní medicína, a.s.	Karasek 1767/1	621 00 Brno CZECH REPUBLIC	Phone: +420-549-124-185 Fax: +420-549-211-460	E-mail: info@biovendor.com sales@biovendor.com Web: www.biovendor.com
AUSTRIA: BioVendor GesmbH	Nußdorfer Straße 20/10	1090 Vienna AUSTRIA	Phone: +43-1-89090-25 Fax: +43-1-89090-2515	E-mail: infoAustria@biovendor.com
GERMANY, SWITZERLAND: BioVendor GmbH	Otto-Hahn-Straße 16	34123 Kassel GERMANY	Phone: +49-6221-433-9100 Fax: +49-6221-433-9111	E-mail: infoEU@biovendor.com
USA, CANADA AND MEXICO: BioVendor LLC	128 Bingham Rd. Suite 1300	Asheville, NC 28806 USA	Phone: +1-828-575-9250 +1-800-404-7807 Fax: +1-828-575-9251	E-mail: infoUSA@biovendor.com